


# Acid-Base Lab

## pH



**Objective:** Students will test a minimum of 15 different substances from “around the house” to determine their pH and then calculate the concentration of  $\text{H}_3\text{O}^{+1}$ .

**Background:** What is pH? pH is a measure of how acidic a substance is and actually stands for “power of hydronium ion.” The pH range is from 0-14 where 0 is acidic and 14 is basic. Many personal hygiene products claim to be “pH balanced.” If they are, what pH “should” they be? Other products you use around the house may be very corrosive (acidic) or caustic (basic). By making a solution of different products, you can easily test the pH. You can then

convert the pH into the concentration of hydronium using the formula,  $\text{pH} = -\log [\text{H}_3\text{O}^{+1}]$ .

### Materials:

- 15 household products to test (can be liquid or solid)
- Watch glasses
- Distilled water
- Stirring rods
- Disposable pipettes
- pH paper

### Procedure:

1. Your watch glasses must be completely CLEAN before use. Rinse THOROUGHLY the watch glass with distilled water. Dry off with paper towel.
2. For solid substances, you will need to make a solution. Place a small amount of the substance onto the watch glass and add just enough distilled water to make a liquid solution. Use the glass stirring rod to mix the water with the solid substances. Proceed to step #4.
3. For liquid substances, use the plastic dropping pipette to place several drops onto a clean watch glass. Proceed to step #4.
4. Take some of the strips of pH paper out of the clear tube. Carefully SEPARATE them into individual strips. Do NOT use “double thick” pH paper!
5. Take ONE piece of pH paper from the vial and tear it in half. You will use one half to test one sample and the other half to test a different sample.
6. Dip one end of the pH paper into the sample. Match its color with the closest shade on the comparison chart. Compare the colors within 30 seconds. Record value on data table (back side of this page).
7. Dispose of the used pH paper in the TRASH.
8. Continue to test samples until done.

## Data Table:

	<b>Sample description</b>	<b>Observed pH value</b>
1.		
2.		
3.		
4.		
5.		
6.		
7.		
8.		
9.		
10.		
11.		
12.		
13.		
14.		
15.		

